

---

# PyCellID

*Release 0.1*

**Clemente, Jose**

**Dec 22, 2021**



## CONTENTS:

<b>1</b>	<b>Installation</b>	<b>1</b>
1.1	Development . . . . .	1
1.2	Initialization . . . . .	1
<b>2</b>	<b>Quick start</b>	<b>3</b>
2.1	Project Structure . . . . .	3
2.1.1	Resource helpers . . . . .	3
2.2	Load your data . . . . .	4
2.2.1	Load your data frame . . . . .	4
2.2.2	From Cell-ID output . . . . .	4
2.3	Inspect your data . . . . .	5
2.3.1	Using Cell-ID features . . . . .	6
2.4	Cell-ID together with PyCellID provide 5 categories of calculated variables: . . . . .	6
2.4.1	View your images . . . . .	8
2.4.2	Use CellsPlotter to inspect images . . . . .	9
<b>3</b>	<b>Basic examples</b>	<b>11</b>
3.1	Prerequisites . . . . .	11
3.2	CellData object instantiation . . . . .	11
3.3	Inspection of data structure . . . . .	12
3.4	Filtering data . . . . .	13
3.5	CellsPlotter . . . . .	14
<b>4</b>	<b>License</b>	<b>19</b>
<b>5</b>	<b>PyCellID Changelog</b>	<b>21</b>
<b>6</b>	<b>pycellid</b>	<b>23</b>
6.1	pycellid package . . . . .	23
6.1.1	Submodules . . . . .	23
6.1.2	pycellid.core module . . . . .	23
6.1.3	pycellid.images module . . . . .	25
6.1.4	pycellid.io module . . . . .	26
6.1.5	Module contents . . . . .	27
<b>7</b>	<b>Indices and tables</b>	<b>29</b>
	<b>Python Module Index</b>	<b>31</b>
	<b>Index</b>	<b>33</b>



## INSTALLATION

PyCellID can be installed using pip from PyPI.

```
$ pip install pycellid
```

### 1.1 Development

We welcome all kinds of contributions. You can build the development environment locally with the following commands:

```
$ git clone git@github.com:pyCellID/pyCellID.git
$ cd pyCellID
$ python3 -m venv venv
$ source venv/bin/activate
$ pip install -e .
$ pip install -r requirements/dev.txt
```

Run the tests with pytest:

```
$ pytest -v tests/
```

Or run the full checks with tox:

```
$ tox -r
```

### 1.2 Initialization

```
>>> import pycellid
>>> pycellid.__version__
'0.0.18'
>>> # Build or load your data, inspect your images and make plots.
>>> from pycellid.core import CellData, CellsPlotter
>>> import pycellid.io as ld # Build or load data frame.
>>> from pycellid import images # Get a 2-D array representing your images.
```



## QUICK START

### 2.1 Project Structure

PyCellID is an Open Source python package designed for recursively navigating a path and tracking down data tables and metadata (mapping), returning a unique CellData object. It relates easily the data with the associated images.

The power of PycellID will be maximized when working with images of time series.

file structure:

```
~folder_location
|
|--- my_experiment/
|   |
|   |---Position01/
|   |   |---mapping.txt
|   |   |---out_all.txt
|   |   |
|   |---Position02/
|   |   |--- ...
|   |   |
|   |---Position03/
|   |   |--- ...
|   |   |
|   |---Position32/
|   |   |--- ...
|   |   |
|   |---channel-x_Position01_time01.tif
|   |---channel-x_Position01_time01.tif.out.tif
|   |--- ...
|   |---channel-z_PositionW_timeY.tif
```

#### 2.1.1 Resource helpers

```
[1]: # Build or load your data, inspect your images and make plots.
    from pycellid.core import CellData, CellsPlotter
    # Build or load data frame.
    import pycellid.io as ld
    # Get a 2-D array representing your images.
    from pycellid import images
    import matplotlib.pyplot as plt
```

## 2.2 Load your data

### 2.2.1 Load your data frame

Instantiate a CellData from your data table and the path to your images.

### 2.2.2 From Cell-ID output

```
[2]: df = CellData.from_csv('samples_cellid')
```

```
[3]: display(df)
```

	pos	t_frame	ucid	cellID	time	xpos	ypos	a_tot	num_pix	\									
0	1	0	100000000000	0	0	29	725	527.0	527										
1	1	1	100000000000	0	0	31	725	614.0	614										
2	1	2	100000000000	0	0	31	724	639.0	639										
3	1	3	100000000000	0	0	31	725	688.0	688										
4	1	4	100000000000	0	0	24	721	457.0	457										
...	...	...	...	...	...	...	...	...	...										
18184	3	12	300000000720	720	0	802	736	147.0	147										
18185	3	12	300000000721	721	0	840	707	391.0	391										
18186	3	12	300000000722	722	0	915	827	587.0	587										
18187	3	12	300000000723	723	0	986	336	119.0	119										
18188	3	12	300000000724	724	0	1014	84	120.0	120										
	fft_stat	...	f_nucl_tag6_tfp	f_nucl_tag6_yfp	f_local_bg_cfp	\													
0	0.387982	...	1675334.0	54200.0	378.6183														
1	0.443545	...	1473843.0	54020.0	383.6143														
2	0.500724	...	1451150.0	46463.0	379.6064														
3	0.591332	...	1566434.0	47871.0	390.2697														
4	0.085026	...	1204067.0	52761.0	372.6058														
...	...	...	...	...	...														
18184	0.556508	...	1237009.0	39586.0	520.8205														
18185	0.040324	...	1053566.0	44035.0	477.8000														
18186	0.672861	...	1534440.0	58636.0	508.6136														
18187	0.494793	...	469596.0	36213.0	530.8571														
18188	0.827304	...	526434.0	127674.0	453.5542														
	f_local_bg_rfp	f_local_bg_tfp	f_local_bg_yfp	f_local2_bg_cfp	\														
0	240.4300	12571.83	312.5942	372.9697															
1	240.4676	12212.17	319.5286	377.4458															
2	241.2868	11988.62	322.5431	376.9560															
3	436.6727	11646.41	330.3378	391.3493															
4	241.3909	11593.46	330.6237	367.9105															
...	...	...	...	...															
18184	254.4151	10049.78	612.9800	522.1818															
18185	248.6977	11320.68	576.1446	463.5887															
18186	250.6716	9132.60	578.7262	509.0000															
18187	250.7143	11290.47	423.6049	536.5789															
18188	249.6765	10899.60	716.8197	460.7500															
	f_local2_bg_rfp	f_local2_bg_tfp	f_local2_bg_yfp																
0	241.2194	12523.050	310.0760																
1	240.1235	12138.300	317.3976																
2	242.0784	11993.090	323.3418																

(continues on next page)



(continued from previous page)

```

3          437.9252          11549.440          331.3429
4          242.5352          11503.960          329.0977
...
18184       254.9767          10533.790          621.8302
18185       250.3025          11368.400          571.9272
18186       250.9851           9129.724          572.0000
18187       251.0638          11238.870          423.4595
18188       248.6575          11796.740          903.9592

```

```
[18189 rows x 125 columns]
```

## 2.3 Inspect your data

The development team decided to use `pandas` library as backend because of its syntax and its extensive documentation. The idea is to make you feel you are working with a pandas object, but with the flexibility of having access to your experimental images.

*You would be able to choose from different backends in future versions.*

```
[4]: df.describe()
```

```

[4]:
count      pos      t_frame      ucid      cellID      time  \
count  18189.000000  18189.000000  1.818900e+04  18189.000000  18189.0
mean      1.966078      6.216944  1.966078e+11      317.438947      0.0
std       0.793782      3.726506  7.937819e+10      204.508054      0.0
min       1.000000      0.000000  1.000000e+11      0.000000      0.0
25%       1.000000      3.000000  1.000000e+11      150.000000      0.0
50%       2.000000      6.000000  2.000000e+11      297.000000      0.0
75%       3.000000      9.000000  3.000000e+11      452.000000      0.0
max       3.000000     12.000000  3.000000e+11      966.000000      0.0

count      xpos      ypos      a_tot      num_pix      fft_stat  \
count  18189.000000  18189.000000  18189.000000  18189.000000  18189.000000
mean      683.000935      529.324097      441.458959      441.460718      0.265314
std      396.587704      289.357639      241.039776      241.038580      0.225960
min       9.000000      9.000000      98.000000      98.000000      0.018116
25%      327.000000      285.000000      277.000000      277.000000      0.068958
50%      681.000000      538.000000      384.000000      384.000000      0.202481
75%     1025.000000      779.000000      555.000000      555.000000      0.420594
max     1384.000000     1032.000000     1499.000000     1499.000000      2.200194

count      ...  f_nucl_tag6_tfp  f_nucl_tag6_yfp  f_local_bg_cfp  f_local_bg_rfp  \
count      ...      1.818900e+04      1.818900e+04      18189.000000      18189.000000
mean      ...      1.258965e+06      1.087345e+05      462.553463      258.523334
std      ...      4.217068e+05      1.837672e+05      42.853657      50.494524
min      ...      2.100890e+05      7.107000e+03      314.227900      232.300000
25%      ...      9.154950e+05      2.818300e+04      439.835400      244.754400
50%      ...      1.177204e+06      4.325000e+04      463.130400      248.077800
75%      ...      1.592472e+06      1.144470e+05      486.517200      251.027400
max      ...      2.342580e+06      2.358444e+06      734.365400      615.073200

count      f_local_bg_tfp  f_local_bg_yfp  f_local2_bg_cfp  f_local2_bg_rfp  \
count      18189.000000      18189.000000      18189.000000      18189.000000
mean      11753.645404      502.869092      459.888134      258.526309
std      1378.053509      349.180950      44.225006      50.631219

```

(continues on next page)

(continued from previous page)

min	6984.529000	0.000000	316.946400	231.608700
25%	10809.180000	350.870100	435.564800	244.753600
50%	11783.600000	397.611100	458.925400	248.076900
75%	12781.720000	511.681500	484.095600	250.980100
max	15438.730000	6654.856000	971.200000	618.813700

	f_local2_bg_tfp	f_local2_bg_yfp
count	18189.000000	18189.000000
mean	11921.632917	489.371202
std	1297.764306	308.290503
min	0.000000	0.000000
25%	11063.320000	351.651900
50%	11939.820000	397.530600
75%	12887.200000	501.397500
max	15294.890000	7044.286000

[8 rows x 125 columns]

### 2.3.1 Using Cell-ID features

## 2.4 Cell-ID together with PyCellID provide 5 categories of calculated variables:

### 1. General measurments.

```
pos,
cellID,
ucid,
t_frame,
time,
xpos,
ypos,
f_tot,
a_tot,
fft_stat,
perim,
maj_axis,
min_axis,
flag,
rot_vol,
con_vol,
a_vacuole,
f_bg
```

### 2. To calculate membrane proximal fluorecence (for relocalization experiment).

```
f_tot_p1_channels,
a_tot_p1,
f_tot_m1_channels,
a_tot_m1,
f_tot_m2_channels,
a_tot_m2,
```

(continues on next page)

(continued from previous page)

```
f_tot_m3_channels,
a_tot_m3
```

### 3. Information obtained from “nuclear image” type (Variables containing the area and fluorescence of concentric disks of user-defined radius).

```
f_nucl_channels,
f_nucl1_channels to f_nucl6_channels,
a_nucl1 to a_nucl6,
f_nucl_tag1_channels to f_nucl_tag6_channels
```

### 4. More background information.

```
f_local_bg_channels,
a_local_bg,
a_local,
f_local2_bg_channels,
a_local2_bg,
a_local2
```

### 5. More volume measurments.


```
a_surf, sphere_vol,
```

*the final tag \_channels indicates that the variable will be repeated for each illumination type.*

For detailed information, reading [Cell-ID](#) documentation is recommended.

```
[5]: df.plot()
```

```
[5]: <AxesSubplot:>
```



## 2.4.1 View your images

Obtain a numpy array representation of an image. Make a crop, operate o simply plot it.

```
[6]: img = plt.imread("samples_cellid/YFP_Position01_time01.tif")


array = images.box_img(im=img, x_pos=640, y_pos=560, radius=30, mark_center=False)

array

[6]: array([[371., 348., 330., ..., 0., 0., 0.],
           [368., 351., 368., ..., 0., 0., 0.],
           [370., 363., 351., ..., 0., 0., 0.],
           ...,
           [ 0., 0., 0., ..., 0., 0., 0.],
           [ 0., 0., 0., ..., 0., 0., 0.],
           [ 0., 0., 0., ..., 0., 0., 0.]])
```

```
[7]: plt.imshow(array)
```


```
[7]: <matplotlib.image.AxesImage at 0xb50faf0>
```



You can use PycellID accessor to inspect images. + Use data from your dataframe for finding what you are looking for.

```
[8]: # Specify your values.
df.plot(array_img_kws={"channel": "tfp", "n": 50, "criteria": {"a_tot": [0, 500]}})

[8]: <AxesSubplot:>
```




### 2.4.2 Use CellsPloter to inspect images

```
[9]: cells = CellsPloter(df)
```

```
[10]: cells.cells_image()
```

```
[10]: <AxesSubplot:>
```





## BASIC EXAMPLES

### 3.1 Prerequisites

1. Clone the repository running the following command:

```
$ git clone https://github.com/pyCellID/pyCellID
$ cd pyCellID
```

2. create a jupyter notebook called 'example.ipynb'

The data used in the following example can be downloaded from the 'samples\_cellid' folder of the [PycellID](#) repository.

```
[1]: import os
import matplotlib.pyplot as plt
import pandas as pd
import pycellid as ld
```

### 3.2 CellData object instantiation

A CellData object can be instantiated by declaring the parameters:

- path: directory that contains the images segmented by Cell-ID.
- df: a Dataframe containing the variables measured for each cell.

```
[2]: base_dir = os.path.join(".", "samples_cellid")
data_dir = os.path.join(base_dir, "pydata", "df.csv")
data = pd.read_csv(data_dir)
```

```
[3]: cells = ld.CellData(path=base_dir, df=data)
```

## 3.3 Inspection of data structure

Information about the way the data is structured can be obtained through intuitive commands.

- Number of rows times number of columns of the table containing Cell-ID data

```
[4]: cells.size
```

```
[4]: 2255436
```

- Number of rows and columns of the table

```
[5]: cells.shape
```

```
[5]: (18189, 124)
```

- List of the first 15 parameters of data table

```
[6]: cells.columns[0:15].to_list()
```

```
[6]: ['pos',  
      't_frame',  
      'ucid',  
      'cellID',  
      'time',  
      'xpos',  
      'ypos',  
      'a_tot',  
      'num_pix',  
      'fft_stat',  
      'perim',  
      'maj_axis',  
      'min_axis',  
      'a_nucl',  
      'rot_vol']
```

- Description of the statistics of parameters

```
[7]: cells["maj_axis"].describe()
```

```
[7]: count      18189.000000  
mean         28.362692  
std          10.433144  
min           3.343855  
25%          20.955810  
50%          25.280250  
75%          34.729460  
max          92.609380  
Name: maj_axis, dtype: float64
```




## 3.4 Filtering data

You would want to look at your images to rule out out-of-focus cells, find dead cells, or remove those cells that biology tells you to be outliers.

It is possible to present exploratory plots with randomly chosen cells satisfying certain criteria.

```
[8]: # Exploratory plot presenting n = 49 randomly chosen cells
# with an area smaller than 500 pixels.
# Out of focus cells, dead cells, clusters and other
# defects can be seen among the images displayed.
cells.plot(array_img_kws={"channel":"tfp", "n":49, "criteria":{"a_tot":[0, 500]})
```


```
[8]: <AxesSubplot:>
```



Data can also be presented using histograms and other types of plots.

```
[9]: # Histogram of cells in the experiment based on their total area.
cells.a_tot.hist(bins=50, alpha=0.5, legend=True)
```

```
[9]: <AxesSubplot:>
```




Data can be filtered applying thresholds to the measured parameters. In this case we want to differentiate cells with areas greater and smaller than 200 pixels and show examples of them in two plots.

```
[10]: # initialize your figure
fig, axs = plt.subplots(1, 2)
fig.set_size_inches(10, 5)

# Filtering your data
cells[cells["a_tot"]<200].plot(ax=axs[0])
cells[cells["a_tot"]>=200].plot(ax=axs[1])

# titles
plt.suptitle('16 random cells')
axs[0].set_title('a_tot < 200 px')
axs[1].set_title('a_tot >= 200 px')

# customize your output
fig.tight_layout()
```



### 3.5 CellsPloter

PyCellID owns the CellsPloter accessor, in charge of rendering the images. You can use the axes provided by CellPloter with the library of your choice.


CellsPloter defines two useful methods: + cells\_image + cimage

```
[11]: # Defining the accessor
ploter = ld.CellsPloter(cells)
```

cells\_image method allows the user to plot a custom number of cells.

```
[12]: ploter.cells_image({"n":50})
```


```
[12]: <AxesSubplot:>
```



It also let the user to define filtering criteria.


```
[13]: ploter.cells_image({"n":20,"criteria":{"maj_axis":[30,40]}})
```

```
[13]: <AxesSubplot:>
```



```
[14]: ploter.cells_image({"n":20,"criteria":{"maj_axis":[10,20]}})
```

```
[14]: <AxesSubplot:>
```



Data filtering performed in the previous section can be also carried out using method *cells\_image* from CellsPloter.


```
[15]: # initialize your figure
fig, axs = plt.subplots(1, 2)
fig.set_size_inches(10, 5)

a_max = cells['a_tot'].max()

# Filtering your data
ploter.cells_image(array_img_kws={'n':16,'criteria':{'a_tot':[0.0,150]}} ,ax=axs[0])
ploter.cells_image(array_img_kws={'n':16,'criteria':{'a_tot':[150,a_max]}} ,ax=axs[1])

# titles
plt.suptitle('16 random cells')
axs[0].set_title('a_tot < 200 px')
axs[1].set_title('a_tot > 200 px')


# customize your output
fig.tight_layout()
```



On the other hand, *cimage* method is used to plot a single cell in a fluorescence channel based on its unique cell identifier (ucid) and the time stamp (t\_frame) corresponding to the instant of interest.

```
[16]: indentifier = { "channel":"BF", "ucid":100000000045, "t_frame":5}
      box_img_kws = {"radius":30}
      ploter.cimage( indentifier,box_img_kws = box_img_kws)
```

```
[16]: <AxesSubplot:>
```



It can be used to show the time evolution of an individual cell. As an example, a single cell is presented at different times in TFP and CFP channels.

```
[17]: ucid = 100000000047
      box_img_kws = {"radius":30}

      # initialize your figure
      fig, axs = plt.subplots(2, 6)
      fig.set_size_inches(24, 8)
```

(continues on next page)


(continued from previous page)

```

channels = ["TFP", "CFP"]
cmaps = ["Greys", "Blues"]

for j, c in enumerate(channels):
    imshow_kws = None
    imshow_kws = {"cmap": cmaps[j]}
    for i in range(6):
        identifier = { "channel":c, "ucid":ucid, "t_frame":i}
        ploter.cimage( identifier, box_img_kws = box_img_kws, imshow_kws=imshow_kws,
→ax=axis[j,i])
        titles = f'{c}: t_frame = {i}'
        axis[j,i].set_title(titles)

```



## **LICENSE**

### MIT License

Copyright (c) 2021 The PyCellID Python Package

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the “Software”), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED “AS IS”, WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.





**PYCELLID CHANGELOG**



## PYCELLID

## 6.1 pycellid package

### 6.1.1 Submodules

#### 6.1.2 pycellid.core module

Merge and analyze tables of characteristics of cells and find images.

**class** pycellid.core.CellData(path, df)

Bases: object

Collapse your data into a single data frame.

Recursively inspect the path, create a unique identifier per cell, and inspect related images.

##### Parameters

- **\_path** (str) – global path to output CellID tables.
- **\_df** (pandas.DataFrame) – Dataframe (output of CellID) containing all the measured parameters of each cell.

##### Returns

- An instance of CellData containing all the information of
- *microscopy experiment*.

##### Examples

```
>>> from pycellid.core import CellData
>>> df = CellData(
    path = '../my_experiment',
    df = my_dataframe
)
```

**classmethod** from\_csv(path, \*\*kwargs)

Build a data frame from csv files contained in path.

A CellData class will be instantiated.

**get\_dataframe**()

Return a copy of the internal \_df.

**property plot**

Represent set of `cells_image` or numerical data.

For `cimage` method you must specify an identifier `id={}`.

**class** `pycellid.core.CellsPlotter` (*cells*)

Bases: `object`

Accessor to plotter class.

Create a representation of each cell within a grid, inspect an entire image or create a snippet of a single cell. Provide a wrapper of pandas methods for plotting. Return axes to plot.

**cells**

An instance of `CellData` containing all the information of microscopy experiment.

Type *CellData*

**cells\_image** (*array\_img\_kws=None, imshow\_kws=None, ax=None*)

Display a random selection of cells on a square grid.

By default it represents a  $4 \times 4$  matrix chosen at random.

**Returns**

**Return type** `ax` to plot or figure.

**Other Parameters**

- **array\_img\_kws** (*dict*) – Set the `pycellid.images.img_array` parameters.  
`n` : number of cells.  
`channels` : “TFP” or another that you have encoded.
- **imshow\_kws** (*dict*) – If you use matplotlib set equal to `plt.imshow`.
- **ax** – Use your axes to plot.

**cimage** (*identifier, box\_img\_kws=None, imshow\_kws=None, ax=None*)

Show a single cell or complete image.

‘`identifier`’ param is required. Reference to a valid image or position. By default, an image with a size of (1392 X 1040) px will be rendered.

**Parameters identifier** (*path or dict*) – path to an image file `dict = {`  
`“channel”:str, “UCID”:int, t_frame”:int }`.

**Returns**

**Return type** `ax` to plot or figure.

**Other Parameters**

- **box\_img\_kws** (*dict*) – Set the `pycellid.images.box_img` parameters.  
**im** [numpy.array] A full fluorescence microscopy image.  
**x\_pos** [int] *x*-coordinate of the center of the cell of interest.  
**y\_pos** [int] *y*-coordinate of the center of the cell of interest.  
**radius** [int] length (in px) between the center of the image and each edge. Default = 90.  
**mark\_center** [bool] mark a black point. Default = False.
- **imshow\_kws** (*dict*) – If you use matplotlib set equal to `plt.imshow`.

- **ax** – Use your axes to plot.

### 6.1.3 pycellid.images module

Images for PyCellID.

Attention! This module will not provide images. This module provides matrix representations for your analysis or to work with your favorite framework.

`pycellid.images.array_img(data, path, channel='BF', n=16, criteria=None)`

Create a grid of images containing cells which satisfy given criteria.

Resulting image has 'n' instances ordered in a grid of shape 'shape'. Each instance corresponds to a image centered in a cell satisfying provided criteria.

#### Parameters

- **data** (`pandas.DataFrame`) – Dataframe (output of `CellID`) containing all the measured parameters of each cell.
- **path** (`str`) – Path to the directory containing the images associated to `data`.
- **channel** (`str`) – Fluorescence channel of the image. The allowed values are 'BF', 'CFP', 'RFP' or 'YFP'.
- **n** (`int`) – Number of instances composing the grid.
- **criteria** (`dict`) – Dictionary containing the criteria of selection of cells.

**Returns** `iarray` – A grid of n images of cells satisfying given criteria.

**Return type** `numpy.array`

**Raises** `ValueError` – If the number of cells satisfying the selection criteria is less than the number of cells to be shown.

`pycellid.images.box_img(im, x_pos, y_pos, radius=90, mark_center=False)`

Create a single image containing an individualised cell.

The resulting image posses a mark in the center of the individualised cell and a pair of delimiters in the right and bottom edges.

#### Parameters

- **im** (`numpy.array`) – A full fluorescence microscopy image.
- **x\_pos** (`int`) – x-coordinate of the center of the cell of interest.
- **y\_pos** (`int`) – y-coordinate of the center of the cell of interest.
- **radius** (`int`) – lenght (in pixels) between the center of the image and each edge.
- **mark\_center** (`bool`) – mark a black point, default = `False`.

**Returns** `iarray` – An array (image) containing an individualised, center-pinned, cell.

**Return type** `numpy.array`

`pycellid.images.img_name(path, ucid, channel, t_frame=None, fmt='.tif.out.tif')`

Construct the name of an image according to output format of `CellID`.

The returned string contains the path and name of the image.

#### Parameters

- **path** (`str`) – Path to the directory containing images asociated to 'data'.

- **ucid** (*int*) – Unique tracking number.
- **channel** (*str*) – Fluorescence channel of the image. The allowed values are 'BF', 'CFP', 'RFP' or 'YFP'.
- **t\_frame** (*int*) – Time tag of the image.
- **fmt** (*str*) – Format of the image to be read.

**Returns** **name** – Name and path of an image according to the output format of CellID.

**Return type** *str*

### 6.1.4 pycellid.io module

in-out implementations for pyCellID.

`pycellid.io.make_df` (*path\_file*)

Make a dataframe with number tracking 'ucid' and 'position'.

**Parameters** **path\_file** (*str*) – Path to CellID's outall data files.

**Returns** **df** – Dataframe with df['ucid'] unique cell identifier.

**Return type** *pandas.DataFrame*

`pycellid.io.merge_tables` (*path, n\_data='out\_all', n\_mdata='\*mapping'*)

Concatenate tables in the path with pandas method.

Transforms the identifying index of each cell from each data table into a temporal index UCID (Unique Cell Identifier) Disaggregate the columns of morphological measurements into columns by fluorescence channel. It uses the mapping present in the metadata file (mapping).

**Parameters**

- **path** (*str*) – Global path to output 'cellID' tables.
- **n\_data** (*str*) – File name to find each data table.
- **n\_mdata** (*str*) – File name to find metadata tables or mapping\_tags.

**Returns** **df** – Dataframe containing 'cellID' data.

**Return type** *pandas.DataFrame*

### Examples

```
>>> import pycellid.io as ld
>>> df=ld.cellid_table(
    path = '../my_experiment',
    n_data = 'out_all',
    n_mdata = 'mapping'
)
```

`pycellid.io.read_df` (*path\_file*)

Read files with data of fluorescence microscopy experiments.

Create a dataframe with the data and rewrite headers format.

**Parameters** **path\_file** (*str*) – Path to files to be read.

**Returns** **df** – Dataframe with data of fluorescence microscopy experiments.

**Return type** pandas.DataFrame

### 6.1.5 Module contents

pyCellID.

An extension that analyze Cell-ID single-cell.





## INDICES AND TABLES

- `genindex`
- `modindex`
- `search`



## PYTHON MODULE INDEX

### p

- `pycellid`, [27](#)
- `pycellid.core`, [23](#)
- `pycellid.images`, [25](#)
- `pycellid.io`, [26](#)



## INDEX

### A

`array_img()` (in module `pycellid.images`), 25

### B

`box_img()` (in module `pycellid.images`), 25

### C

`CellData` (class in `pycellid.core`), 23

`cells` (`pycellid.core.CellsPloter` attribute), 24

`cells_image()` (`pycellid.core.CellsPloter` method), 24

`CellsPloter` (class in `pycellid.core`), 24

`cimage()` (`pycellid.core.CellsPloter` method), 24

### F

`from_csv()` (`pycellid.core.CellData` class method), 23

### G

`get_dataframe()` (`pycellid.core.CellData` method), 23

### I

`img_name()` (in module `pycellid.images`), 25

### M

`make_df()` (in module `pycellid.io`), 26

`merge_tables()` (in module `pycellid.io`), 26

module

`pycellid`, 27

`pycellid.core`, 23

`pycellid.images`, 25

`pycellid.io`, 26

### P

`plot()` (`pycellid.core.CellData` property), 23

`pycellid`

    module, 27

`pycellid.core`

    module, 23

`pycellid.images`

    module, 25

`pycellid.io`

    module, 26

### R

`read_df()` (in module `pycellid.io`), 26